INSTRUCTIONS FOR USE

The Kapp Cosgrove[®] Valve Retractor

Manufactured by:

Kapp Surgical Instrument, Inc.

4919 Warrensville Center Rd.

Cleveland, Ohio 44128, USA

Telephone: 1-216-587-4400 or 1-800-282-KAPP

Fax: 1-216-587-0411

www.kappsurgical.com e-mail: info@kappsurgical.com

Rx only

Cosgrove is a registered trademark of Kapp Surgical Instrument, Inc.

Cat # KS7000 MAN Rev D (1-21)

Consistent and excellent exposure of the right and left atrium, tricuspid, and mitral valves is fundamental to intracardiac surgical procedure. By using the Kapp Cosgrove® Valve Retractor, the surgeon can gain the optimum exposure necessary for a successful surgical outcome.

DESCRIPTION / INTENDED USE

The Kapp Cosgrove® Valve Retractor provides exceptional right and left atrium, tricuspid and mitral valve exposure while reducing the need for surgical assistance. The Kapp Cosgrove® Valve Retractor requires placement only during the intracardiac portion of the procedure. It is available either mounted on a Cooley Sternum Retractor, large or small, or custom mounted on your existing chest spreader. Recommendations for cleaning and sterilization are included in this manual.

The device is supplied non-sterile and must be cleaned and sterilized before use.

PRE-SURGICAL ASSEMBLY

Position your modified chest spreader as shown (**Figure 1**). Slide the main bracket (**Figure 3**) onto the chest spreader with locking screws to the outside of the chest cavity opening by sliding it onto one of the pre-attached mounting bars (**Figure 2**). Secure the main bracket by tightening the left and right non-removable locking screws (A) (**Figure 3**).

Slide the retractor blade (D) through the universal clamp (C) (**Figure 4**).

Remove the transverse nut (E), then slide universal clamp(s) with attached retractor blade(s) (**Figure 5**) onto the support bar. Rotate the support bar down into the groove in the side support bracket and tighten the transverse nut to secure (**Figure 6**).

Adjust the retractor blades on the main bracket according to the surgeon's preference. Once they are in place, secure the retractor blades by tightening the "wing" nut on the universal clamp(s).

A third retractor blade can be mounted by attaching the side arm bar (F) as follows: Loosen the side arm locking screw (non-removable) on the main bracket (B) (**Figure 3**). Slide the flattened end of the side arm bar between the support bracket and the head of the locking screw. Tighten the locking screw. Slide a universal clamp with attached retractor blade in the bottom opening onto the side arm bar and secure it by tightening the wing nut on the universal clamp (**Figure 6**).

The Kapp Cosgrove® Valve Retractor assembly as shown on the front cover is now complete.

NOTE: All screws on the Kapp Cosgrove® Valve Retractor are non-removable except for the large transverse nut (E) (**Figure 3**).

NOTE: The retractor blades may be easily readjusted during surgery by loosening the wing nut on the universal clamp(s).

CLEANING

Immediately after surgery, rinse the instruments under warm, not hot, water to remove all blood, body fluids and tissue. If delays prevent immediate rinsing, soak the instrument in a solution of water and neutral (pH 7) detergent prior to cleaning.

ULTRASONIC CLEANING

Separate the components and place them into the ultrasonic cleaner. Clean per instrument protocol. Rinse the instruments with water immediately after ultrasonic cleaning.

AUTOMATIC WASHER STERILIZERS

Follow manufacturer's recommendations ensuring instruments are lubricated after last rinse cycle and before sterilization cycle.

According to the Robert Koch Recommendations from September 2012, mechanical cleaning and disinfection should be conducted at temperatures in excess of 93° C.

AUTOClave STERILIZATION

Prior to autoclave sterilization, lubricate all instrument parts which have metal-to-metal contact using surgical lubricants such as instrument milk. Do not use industrial lubricants. Instruments must be unlocked before cleaning to prevent any cracks that may occur from heat expansion during the autoclave cycle. Do not overload the autoclave chamber. The manufacturer's recommended sterilization protocol is:

Sterilizer Type:	Pre-vacuum
Preconditioning Pulses:	3
Minimum Temperature:	132°C
Full Cycle Time:	4 minutes
Minimum Dry Time:	30 minutes
Article Configuration:	Instruments placed in a general surgical case with lid 11 in x 26 in x 3.5 in (28 cm x 66 cm x 9 cm) without brackets or pin mat with the tray wrapped in two layers of 1 ply polypropylene wrap using sequential envelope techniques.

The following cycle should be used when there is an immediate need for the device and there are no other alternatives:

Sterilizer Type:	Pre-vacuum
Preconditioning Pulses:	3
Minimum Temperature:	132°C
Full Cycle Time:	3 minutes
Minimum Dry Time:	N/A
Article Configuration:	Unwrapped
Exposure time:	Minimum 3 minutes

Kapp Surgical Instrument, Inc. confirms that the processing recommendations related to machine cleaning and sterilization by Robert Koch Recommendations from September 2012 may be used. Specifically, Kapp Surgical Instrument, Inc. confirms that sterilization times and temperatures may be extended to 5 minutes and 134° C.

MAINTENANCE REQUIREMENTS

All Kapp Cosgrove® Valve Retractors and attachments should be serviced at least annually. Your retractor should be refurbished by the original manufacturer and original patent holder, Kapp Surgical Instrument, Inc. There is no limit to the number of reuses of this device as long as it is appropriately maintained and serviced by Kapp Surgical Instrument, Inc.

NOTE: All screws on the Kapp Cosgrove® Valve Retractor are non-removable. These screws are peened over to prevent the screws from falling into the surgical field. The only item that is removable is the large transverse nut (E) located on the main bracket (**Figure 3 & 5**).

DO NOT SEND THIS INSTRUMENT OUT TO BE SERVICED TO YOUR LOCAL SURGICAL INSTRUMENT REPAIR COMPANY. PLEASE SEND ALL REPAIRS AND REFURBISHING BACK TO THE ORIGINAL MANUFACTURER:

Kapp Surgical Instrument, Inc.
 4919 Warrensville Center Rd.
 Cleveland, Ohio 44128, USA
 1-800-282-5277 or 1-216-587-4400

INSTRUMENT COMPONENTS

The Kapp Cosgrove® Valve Retractor is made from instrument-grade quality stainless steel. NOTE: This product is normally intended to be used for less than 60 minutes.

WARRANTY

Kapp Surgical instruments are guaranteed to be free from defects in workmanship and material. Any instrument that proves defective in workmanship or material will either be repaired or replaced at our discretion, without charge.

Kapp Surgical Instrument, Inc. warrants the guarantee of this instrument for the life of this product under the normal surgical use for which it was intended.

Illustrations:
Illustraties:
Illustrations:

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