

# Adult Renal Clamps

## Special Cardiovascular & Thoracic Clamps

These CVT Surgical Clamps are ideally suited for clamping of arteries located in the depths and superficial vessels of the surgical field. The double angulation of the clamp permits placement of its handle to one side, allowing an un-obstructed view of the operative field.

The instrument is manufactured in pairs, each being the mirror image of the other, for use to one side or the other, depending upon the circumstance.

The large Adult Renal Clamp (Angled Clamp) is designed to control the common iliac artery, the external or internal iliac arteries during surgery for abdominal aortic aneurysmectomy for aorto-iliac occlusive disease. It can also be used for controlling the distal external iliac artery during aorto-bi-femoral bypass procedures, femoral-popliteal reconstruction, profundaplasties, or femoral embolectomies.

### Adult Renal Clamps (Angled Clamps)

