




**Fazio Colorectal
Cleveland Clinic Intestinal Dilators (Bougies)**



Dr. Victor Fazio, Chairman of the Colorectal Department at the Cleveland Clinic Foundation, originally designed and endorses these dilator / Bougies. The Bougies are made of a surgically approved, latex free, rubber material, with a consistency from flexible to firm, and are hollow. They are **13 inches long** and vary in diameter sizes from **9mm to 32mm**.

A small set of 7 Bougies is available from 9mm to 15mm and a large set of 6 Bougies is available from 18mm to 32mm. A complete set of all 13 Bougies is also available. See catalog listing below.

These dilators are useful for sounding and dilating the lumina of various parts of the intestinal tract at the time of operation, and at the time of patient examination in the outpatient department. They are also useful for assessing a true diameter of the contracted bowel for facilitating the surgical reunion of the divided ends of the intestinal tract, both by conventional suture methods and by stapling methods.

They are useful also in assessing sizes of strictured rectal outlet changes, or for overcoming stenosis of such small diameter that uncomfortable symptoms of partial obstruction are produced.

Catalog No.	KS-6000	Complete Set 13 sizes
	KS-6001	Small Set 7 sizes
	KS-6002	Large Set 6 sizes
Small Set:	KS-6003	9mm
	KS-6004	10mm
	KS-6005	11mm
	KS-6006	12mm
	KS-6007	13mm
	KS-6008	14mm
	KS-6009	15mm
Large Set:	KS-6010	18mm
	KS-6011	20mm
	KS-6012	23mm
	KS-6013	26mm
	KS-6014	29mm
	KS-6015	32mm

KAPP SURGICAL INSTRUMENT, INC

Worldwide Sales, Service & Manufacturing

4919 Warrensville Center Road - Cleveland, Ohio 44128

216-587-4400 - Fax: 216-587-0411 - 1-800-282-KAPP

Web Site: www.kappsurgical.com - Email: info@kappsurgical.com