

Kapp Bulldog Clamps

Special Cardiovascular & Thoracic Clamps

Unlike all other Bulldog Clamps, these clamps allow the surgeon complete control of the pressure required as needed. The vast range of pressure (from very weak to very strong) allowed by this Bulldog can be attributed to its unique design and variable spring. These clamps are manufactured straight 1 3/4" and 2 1/4" long, also 2" long with a right angle 1/4" tip if desired.

Kapp Bulldog Clamps (Shown Actual Size)

Cat # KS-3003
Right Angle 2"



Cat # KS-3002
Long 2 1/4"
Straight Jaws



Cat # KS-3001
Short 1 3/4"
Straight Jaws

