



Vasectomy Instruments

These “no scalpel” Vasectomy instruments are designed to simplify the operation by eliminating the scalpel. The forcep, a sharp pointed and curved instrument is designed to puncture the scrotal sac, pierce the vas where it is then secured by the Fixation Ring Clamp. Then the vas is clipped and divided. The I.D. of the eyelet is approximately 3.5mm. These instruments are made of high quality stainless steel and secure lock ratchets with no sharp edges.



FORCEP

The sharp pointed, and curved hemostat is designed to puncture the scrotal sac, pierce the vas, where it is clipped and divided.

Cat # K-9554



FIXATION RING CLAMP

The fixation ring clamp is used to isolate and immobilize the vas within the scrotum prior to clipping.

Cat # K-9555