



Ancillary Universal Stabilizers

Rigid Universal Stabilizer (Malleable)




Rigid Stabilizer Attachment.
Cat. #KS-8300



Rigid Stabilizer mounted on a common chest spreader used for coronary surgery.

Suction & Adjustable Stabilizers



Suction Stabilizer Tube.
To remove fluid and or flush CO₂.
Cat. #KS-8302




Suction Stabilizer Base
with suction thru stabilizer feet.
Cat. #KS-8303




Universal Ball and Socket
Adjustable Stabilizer.
Cat. #KS-8304

Flexible Universal Stabilizer



Actual Flexible Stabilizer.
Cat. #KS-8400
(Used for coronary surgery)



Stabilizer Distal Fork can be customized to the surgeon's specification.



Flexible Stabilizer can be mounted on any surgical chest spreader.