



# Sterile Adhesive Organizer

**FAST ■ EASY ■ SECURE ■ VERSATILE**



**Self adhesive STRIP-T'S™ hold and organize all surgical tubes, scopes and instruments any place on the surgical field.**

- Economical - Easy to use
- Latex Free
- Saves pre-op time
- Versatile, can be used on every specialty surgery
- Several sizes and configurations
- Iridescent color for laser surgery
- Sterile shelf life 2 years


**STRIP-T'S™** is the least expensive apparatus to hold implements and is less expensive than tape, pouches or drapes. The **STRIP-T'S™** reduces the cost of replacing tubes, costly disposable suctions, surgical operative scopes and bovies, which tend to fall from the operative field.



**SURGICAL  
CONCEPTS, INC.**

a division of KAPP SURGICAL INSTRUMENT, INC.

Worldwide Sales, Service & Manufacturing  
4919 Warrensville Center Rd., Cleveland, Ohio 44128  
Phone: 216-587-4400 / Fax: 216-587-0411 / 800-282-5277  
Web site: [www.kappsurgical.com](http://www.kappsurgical.com) / Email: [info@kappsurgical.com](mailto:info@kappsurgical.com)




**FAST ■ EASY ■ SECURE ■ VERSATILE**

**Conveniently Holds and Organizes Surgical Implements In Every Surgical Procedure**



LAPAROSCOPIC CHOLECYSTECTOMY



LAPAROSCOPIC CHOLECYSTECTOMY



GENERAL SURGERY  
Cat. #SC-5000-3Y




ORTHOPEDIC SURGERY  
Arthroscopy & total knee replacement shown with bovie scratch pad.



ORTHOPEDIC SURGERY  
Total hip replacement.

**STRIP-T'S™ STANDARD CONFIGURATIONS**



Cat. #SC-5000-3Y (Yellow) or SC-5000-3R (Red)

For use on all surgical procedures. Three self-closing Hook and Loop compartments plus one open Hook and Loop section. Two Hook and Loop pencil strips are provided to quickly wrap around your instrument or bovie and stick to the open Hook and Loop section. Scratch emery attached. Measures 2" x 12". Each Hook and Loop pencil strip measures 5/8" x 2-1/2".



Cat. #SC-5000-4Y (Yellow)

For general use on all surgeries. Four self-closing Hook and Loop compartments secure cords, tubes, scopes and instruments. Measures 2" x 10".



Cat. #SC-5000-4R (Red)

Note: Custom lengths and designs are available with 2,500 unit minimum order.



**SURGICAL CONCEPTS, INC.**

a division of KAPP SURGICAL INSTRUMENT, INC.

Worldwide Sales, Service & Manufacturing  
4919 Warrensville Center Rd., Cleveland, Ohio 44128  
Phone: 216-587-4400 / Fax: 216-587-0411 / 800-282-5277  
Web site: [www.kappsurgical.com](http://www.kappsurgical.com) / Email: [info@kappsurgical.com](mailto:info@kappsurgical.com)