




**KAPP SURGICAL  
INSTRUMENT, INC.**

***Loop Aorta Clamp***

## ***Loop Aorta Special CVT-Soft Jaw Aorta Clamp***

We have devised a light, strong clamp for use on the distal ascending aorta. Its advantages are that it may be used with "cushioned jaws" or similar padded inserts, and it conforms to the depth of the superior mediastinum. The clamp lies flat against the neck and does not protrude near the operative field.



Dr. Floyd Loop  
CVT Surgeon  
CEO of Cleveland Clinic Foundation

Catalog # KS-7700  
Catalog # KS-7700 XL

### **KAPP SURGICAL INSTRUMENT, INC.**

Worldwide Sales, Service & Manufacturing  
4919 Warrensville Center Road  
Cleveland, Ohio 44128

(216) 587-4400 FAX (216) 587-0411  
1-800-282-KAPP

Web Site: [www.kappsurgical.com](http://www.kappsurgical.com)  
Email: [info@kappsurgical.com](mailto:info@kappsurgical.com)

