


Martinez Clamps

Special Cardiovascular & Thoracic Clamps

These CVT Surgical Clamps are ideally suited for clamping of arteries located in the depths and superficial vessels of the surgical field. The double angulation of the clamp permits placement of its handle to one side, allowing an unobstructed view of the operative field.

The instrument is manufactured in pairs, each being the mirror image of the other, for use on one side or the other, depending on the circumstance. The Martinez Clamp (Angled Clamp) is designed to control the internal, external, and common carotid arteries during carotid endarterectomies, for the control of brachial artery during thrombectomy and/or repair following cardiac catheterization, and for profunda femoris artery control during femoral artery procedures.

Martinez Clamps Angled Clamp with Cooley Teeth



**KAPP SURGICAL
INSTRUMENT, INC.**

4919 Warrensville Center Rd. • Cleveland, Ohio 44128
800-282-KAPP • Phone: 216-587-4400 • Fax: 216-587-0411
www.kappsurgical.com • Email: info@kappsurgical.com

Rev. 11-21 v2