


Rosenkranz™ Pediatric Open Heart Retractor

A unique retractor with attachable, self-retaining atrial and ventricular retractor blades that provide excellent and consistent exposure of all heart chambers during heart valve replacement for congenital heart defect surgery. The Rosenkranz™ Pediatric Open Heart Retractor also enhances visualization of the intra-cardiac anatomy.

FEATURES

- Versatile components are easily positioned or removed from the sternal retractor during any portion of the procedure
- Universally adjustable retractor blades accommodate any anatomic situation for greater visualization while reducing the need for surgical assistance

For excellent and consistent heart valve exposure, the Rosenkranz™ Pediatric Open Heart Retractor sets the standard for pediatric and small child use.



Patented



SYSTEMS

- KS 7500** Rosenkranz™ Pediatric (0-1 year) W/Sternal Retractor
- KS 7501** Rosenkranz™ Small Child (1-5 years) W/Sternal Retractor
- KS 7502** Rosenkranz™ Pediatric attachment mounted on your own Retractor
- KS 7529** Rosenkranz™ Pediatric attachment only

COMPONENTS

- KS 7520** Rosenkranz™ Main Bracket
- KS 7521** Rosenkranz™ Universal Clamp
- KS 7522** Rosenkranz™ Small Retractor Blade
- KS 7523** Rosenkranz™ Deep Retractor Blade
- KS 7524** Rosenkranz™ Wire Basket Retractor
- KS 7525** Rosenkranz™ Side Arm Attachment
- KS 7526 P** Rosenkranz™ Universal Malleable Retractor Blade (copper – K-217) Pediatric 3 1/8"L x 3/8"W
- KS 7526 C** Rosenkranz™ Universal Malleable Retractor Blade (copper – K-216) Child 4 3/4"L x 3/8"W
- KS 7527** Rosenkranz™ Pediatric Sternal Retractor 19mmx 21mmx 75mm
- KS 7528** Rosenkranz™ Small Child Sternal Retractor 21mmx 45mmx 140mm

CUSTOM PARTS and SERVICES

- KS 7503** Dual mounting bars added to your own pediatric or small child sternal retractor to mount Rosenkranz™ Pediatric Open Heart Retractor attachment

